99207 Rev. C 2008-09 • Änderungen vorbehalten • Subject to alterations • Sous réserve de modifications • Con riserva di apportare modifiche • Внесены изменения • B03M0ЖНы изменения

# **Riester**

The familiar way

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**Diagnostische  
Instrumente**

Instructions

**Diagnostic Instruments**

Mode d'emploi

**Instruments**

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**Instrumentos**

**diagnósticos**

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# **1. Important information - read prior to start-up**

You have acquired a valuable RIESTER diagnostic set manufactured in compliance with Directive 93/42/EEC for medical products and subject to continuous stringent quality control, whose excellent quality will ensure reliable diagnoses.

Please read these Operating Instructions carefully prior to start-up and keep in a safe place.

Should you have any queries, please contact the Company or your RIESTER Agent who will be pleased to assist you. For addresses see last page of these Operating Instructions. The address of your authorised RIESTER Agent will be supplied to you on request.

Please note that any instruments described in these Operating Instructions are only suited for application by trained operators.

Please also note that correct and safe operation of instruments will only be guaranteed when RIESTER instruments and accessories are used throughout.

## **2. Battery handles and start-up**

### **2.1. Purpose**

RIESTER battery handles described in these Instructions for Use supply the instrument heads with power (the lamps are included in appropriate instrument heads), also serving as a bracket.

### **2.2. Readiness for operation (insertion and removal of batteries)**

Turn off instrument head from handle in counter-clockwise direction. Insert two commercial type "AA" Mignon alkaline batteries of 1.5 V (IEC standard reference LR6) into the case of the handle with the plus poles towards the upper section of the handle.

#### **Warning:**

- Should the unit not be used for an extended period of time or whilst travelling, remove batteries from handle.
- Insert new batteries when light intensity of the unit is reduced, thus affecting examination.
- For maximum light yield it is recommended to always insert new high-quality batteries on replacement.
- Ensure that no fluid or condensation penetrates into the handle.

## Disposal

Please note that batteries are subject to separate disposal. For details ask your local authority and/or your environmental officer.

### 2.3. Attachment of instrument heads

Turn instrument head in clockwise direction on to the handle. The instrument head may be fixed in the required position by turning the ring on the instrument head up to the stop.

### 2.4. Starting and stopping

ON The clip on the handle is equipped with an On/Off slide. When pushing the slide up, the unit is switched on, when pushing it down, the unit is off.

OFF

### 2.5. Care

#### Cleaning and/or disinfection

Battery handles may be cleaned externally with a moist cloth. In addition, the following disinfectants may be used for disinfection:

Aldehyde (formaldehyde, glutaraldehyde, aldehyde derivatives), surfactants or alcohols.

When using these substances, the manufacturer's instructions must be strictly complied with.

Means for cleaning or disinfection may be a soft, possibly lint-free cloth or Q-tips.

#### Warning

Never immerse handles in fluid.

Ensure that no fluid penetrates into the interior of the enclosure.

#### Sterilisation

According to the ruling teaching (Tübingen Medical Products Test

Centre), sterilisation is only specified for surgery.

The battery handle is not used in surgery. Therefore, cleaning or dis-infecting will be fully adequate.

### Classification

Application type B

### 3. Otoscope and accessories

#### 3.1. Purpose

RIESTER otoscopes described in these Instructions for Use have been produced for lighting and examination of the auditory canal, combined with a RIESTER ear speculum.

#### 3.2. Insertion and removal of ear speculum

Position the selected speculum on the chromium plated metal socket of the otoscope. Turn speculum to the right until a resistance is felt. The size of the speculum is marked on the reverse.

#### 3.3. Introduction of external instruments into the ear

When intending to introduce external instruments into the ear (such as forceps), turn magnifying glass (approx. 2.5X enlargement) on otoscope head in counter-clockwise direction. Replace cover lens in reverse direction.

#### 3.4. Pneumatic test

In order to perform a pneumatic test (= examination of the ear drum), you will require a bulb which is not included in the normal scope of supply but may be ordered separately (see Spare parts and accessories).

Take metal connector which is not included in the normal scope of supply but may be ordered separately (see Spare parts and accessories) and insert in recess provided on the side of the otoscope head. Attach hose of bulb to connector. Carefully introduce the required air volume into the auditory canal.

### 3.5. Replacement of lamp

#### ri-mini®

Remove instrument head from battery handle.

The lamp is in the bottom section of the instrument head. Remove lamp, by using your thumb and forefinger or a suitable tool, from instrument head. Firmly insert new lamp.

#### pen-scope®

Screw off the metal fitting. Release lamp by turning in counter-clockwise direction.

Insert new lamp by screwing in in clockwise direction and screw on the metal fitting again.

### 3.6. Specification of lamp

ri-mini® otoscope	2,5 V average life 15 h	750 mA
pen-scope® otoscope	2,7 V/HL2,5 V average life 10 h	300 mA

### 3.7. Care

#### Cleaning and/or disinfection

The otoscope may be cleaned externally with a moist cloth.

In addition, the following disinfectants may be used for disinfection:

Aldehyde (formaldehyde, glutaraldehyde, aldehyde derivatives), or surfactants.

In addition, all instrument components, with the exception of the magnifying glass, may be disinfected with alcohol.

When using these substances, the manufacturer's instructions must be strictly complied with.

Means for cleaning or disinfection may be a soft, possibly lint-free cloth or Q-tips.

#### Warning

Never immerse the otoscope in fluid.

Ensure that no fluid penetrates into the interior of the enclosure.

#### Sterilisation

##### a) Otoscope

According to the ruling teaching (Tübingen Medical Products Test

Centre), sterilisation is only specified for surgery. Due to the otoscope not being used in surgery, cleaning or disinfecting will be fully adequate.

##### b) Reusable ear specula

Although sterilisation is not required, as specified in a), the ear speculum may nevertheless be sterilised at 134° C for a dwell time of 10 minutes in a steam steriliser.

### 3.8. Spare parts and accessories

#### Reusable ear specula

Item No. 10775	Pack of 10 St.	2 mm, black
Item No. 10776	Pack of 100 St.	2 mm, black
Item No. 10777	Pack of 500 St.	2 mm, black
Item No. 10778	Pack of 1000 St.	2 mm, black
Item No. 10779	Pack of 10 St.	2,5 mm, black
Item No. 10780	Pack of 100 St.	2,5 mm, black
Item No. 10781	Pack of 500 St.	2,5 mm, black
Item No. 10782	Pack of 1000 St.	2,5 mm, black
Item No. 10783	Pack of 10 St.	3 mm, black
Item No. 10784	Pack of 100 St.	3 mm, black
Item No. 10785	Pack of 500 St.	3 mm, black
Item No. 10786	Pack of 1000 St.	3 mm, black
Item No. 10789	Pack of 10 St.	4 mm, black
Item No. 10790	Pack of 100 St.	4 mm, black
Item No. 10791	Pack of 500 St.	4 mm, black
Item No. 10792	Pack of 1000 St.	4 mm, black
Item No. 10795	Pack of 10 St.	5 mm, black
Item No. 10796	Pack of 100 St.	5 mm, black
Item No. 10797	Pack of 500 St.	5 mm, black
Item No. 10798	Pack of 1000 St.	5 mm, black

#### Spare lamps

Item No. 10488	Pack of 6 St. lamps 2,7 V
Item No. 10489	Pack of 6 St. lamps HL 2,5 V
	<b>for pen-scope® Otoscope</b>
Item No. 10600	Pack of 6 St. halogen lamps 2,5 V
	<b>for ri-mini® Otoscope</b>

#### Other spare parts

Item No. 10960	Bulb for pneumatic test
Item No. 10961	Connector for pneumatic test
Item No. 10963	Magnifying glass for ri-mini® and pen-scope® Otoscope, black
Item No. 10962	Magnifying glass for ri-mini® and pen-scope® Otoscope, slate gray
Item No. 10477	Magnifying glass for ri-mini® and pen-scope® Otoscope, blue
Item No. 10964	Magnifying glass for ri-mini® and pen-scope® Otoscope, saffron
Item No. 10965	Magnifying glass for ri-mini® and pen-scope® Otoscope, green

## 4. Ophthalmoscope and accessories

### 4.1. Purpose

RIESER ophthalmoscopes described in these Instructions for Use have been designed for the examination of the eye and its background.

### 4.2. Lens wheel and correcting lenses

The correcting lenses may be adjusted on the lens wheel. The following correcting lenses are available: diopters 0 to + 20 and 0 to - 20. Readings will be displayed on a lit panel. Plus values are displayed in black digits, minus values in red digits.

### 4.3. Apertures and filters


The following apertures and/or filters may be selected by the aperture and filter wheel:

Aperture	Function
● Small circle	
◐ and semi-circle:	For reduction of reflexes of small pupils
● Large circle:	For standard fundus examination
⊙ Fixation star:	For definition of central and eccentric fixation

**The following filters are only available in the ri-mini® ophthalmoscope:**

○ Red-free filter (green filter)	To increase contrast for assessment of changes in fine vessels, i.e. retinal haemorrhages
----------------------------------	---

### 4.4. Lamp replacement

Remove instrument head from battery handle.

The lamp is in the bottom section of the instrument head. Remove lamp, by using your thumb and forefinger or a suitable tool, from the instrument head. Insert new lamp with the pin on the lamp fitted in the recess provided on the head of the case.

### 4.5. Specification of lamp

ri-mini®	Ophthalmoscope 2,5 V average life 15 h	750 mA
pen-scope®	Ophthalmoscope 2,7 V average life 15 h	300 mA

### 4.6. Care

#### Cleaning and/or disinfection

The ophthalmoscope may be cleaned externally with a moist cloth. In addition, the following disinfectants may be used for disinfection: Aldehyde (formaldehyde, glutaraldehyde, aldehyde derivatives), or surfactants. All instrument components with the exception of the lens wheel may be disinfected with alcohol. When using these substances, the manufacturer's instructions must be strictly complied with. Means for cleaning or disinfection may be a soft, possibly lint-free cloth or Q-tips.

### **Warning**

Never immerse the ophthalmoscope in fluid.  
Ensure that no fluid penetrates into the interior of the enclosure.

### **Sterilisation**

According to the ruling teaching (Tübingen Medical Products Test Centre), sterilisation is only specified for surgery.  
Due to the ophthalmoscope not being used in surgery, cleaning or disinfecting will be fully adequate.

## **4.7. Spare parts and accessories**

### **Spare lamps**

Item No. 10473	Pack of 6 standard lamps of 2.7 V <b>for pen-scope® Ophthalmoscope</b>
Item No. 10605	Pack of 6 St. halogen lamps of 2.5 V <b>for ri-mini® Ophthalmoscope</b>

## **5. Maintenance**

These instruments and their accessories do not require any specific maintenance.  
Should an instrument have to be examined for any specific reason whatsoever, please return it to the Company or an authorised RIE-STER dealer in your area. Addresses to be supplied on request.

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## **Riester offers a large selection of products in the areas of**

Blood pressure measuring devices | Instruments for ENT,  
Ophthalmological instruments | Dermatological instruments  
| Thermometers | Stethoscopes | Head mirrors, Head lights,  
Examination lights | Laryngoscopes | Gynaecological instru-  
ments | Percussion hammers | Tuning forks | Products for  
blood stasis | Pulmonary pressure measuring devices |  
Dynamometers | Pressure infusion instruments | Veterinary  
instruments | Doctor's cases and bags

**Detailed descriptions of the products can be found in the respective sections of the omnibus edition catalogue (Order No. 51232-50). Or online under [www.riester.de](http://www.riester.de).**

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## **Rudolf Riester GmbH**